- DE** Bedienungsanleitung
- EN** Operation Manual
- FR** Mode d'emploi
- IT** Istruzione d'uso
- ES** Instrucciones de uso

HighSurg 30

Motor System for Medical Surgery

DE

EN

FR

IT

ES

Vorwort

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Foreword

Congratulations on your purchase of a NOUVAG AG product. Thank you for the confidence shown in our products. Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years.

You will find the conformity statement and list of authorized service representatives attached.

- **Please read instructions carefully before operating!**

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Table of contents

1	Product description	2
1.1	Intended use and operation	2
1.2	Contraindications	2
1.3	Technical data, HighSurg 30	2
1.4	Ambient conditions	2
1.5	Warranty coverage	2
2	Explanation of symbols	3
3	Safety information	3
3.1	EMC Manufacturer's Declaration of Conformity	3
3.2	Integrated peristaltic tube pump	3
3.3	Modification and misuse	4
3.4	Essential requirements	4
3.5	During use	4
4	Scope of delivery	5
5	Device overview	6
6	Startup	7
6.1	Device setup	7
6.2	Connection to the power supply	7
6.3	Device preparation	8
6.4	Assembly of external irrigation system	10
7	Operation	11
7.1	Switching the device on and off	11
7.2	Overview: Control elements on the operation panel	11
7.3	Overview: standard display	12
7.4	Adjusting the programs	13
7.4.1	<i>Step 1: Selecting handpiece or transmission ratio</i>	13
7.4.2	<i>Step 2: Calibrating handpieces</i>	14
7.4.3	<i>Step 3: Setting the speed</i>	15
7.4.4	<i>Step 4: Setting the torque</i>	15
7.4.5	<i>Step 5: Setting the pump supply quantity</i>	16
7.5	Torque limiter	16
7.6	Programme memory	16
7.7	Configuration Menu	17
7.8	Operation using the Vario pedal	21
7.9	Functional check	22
7.9.1	<i>Electronic motor</i>	22
7.9.2	<i>Pump</i>	22
7.9.3	<i>Rotational direction of the electronic motor</i>	22
7.9.4	<i>Program</i>	22
8	Cleaning, disinfection and sterilization	23
8.1	Control unit and pedal	23
8.2	OTO-Drill handpiece with integrated Electronic Motor 16	23
8.3	Electronic motor 21	24
8.4	Tubing set, Ref. 6024 and Ref. 6025	25
8.5	Handpiece cradle	25
9	Maintenance	26
9.1	Replacing the control unit fuse	26
9.2	Safety inspections	26
10	Malfunctions and troubleshooting	27
11	Spare parts list with order numbers	30
12	Information on disposal	30

EN

1 Product description

1.1 Intended use and operation

The HighSurg 30 in connection with an electronic motor and a corresponding handpiece (separate medical product) is applied in general surgery (drilling and milling) by 1 : 1 handpieces with FG-shank, as well as in further different medical fields:

- Plastic Surgery: by K-Wire Handpiece, Dermatomes, Tattooing handpiece
- Spinal Surgery: by Spinal Burr Handpieces, Saws, Multiplier Handpiece.
- Head and Neck Surgery/Skull Surgery/ENT Surgery: by straight and angled Drilling Handpieces 1 : 1, Saws, Craniotomes, Perforators as well as by straight and angled Shaver Handpieces (Micro Debridors).
- Orthopedics/Arthroscopy: by straight and angled Drilling Handpieces 1 : 1 and 4 : 1, K-Wire Handpiece, Saws as well as by straight and angled Shaver Handpieces (Micro Debridors).

The patient group is not limited by age, weight or gender.

The setting up and operation of the HighSurg 30 must be carried out by the physician or by a qualified especially trained medical personnel.

1.2 Contraindications

Relative or absolute contraindications can arise from the general medical diagnose, or in special cases by a significantly increased risk to the patient through the use of motor-driven systems.

Relevant cases in the literature must be taken into consideration.

1.3 Technical data, HighSurg 30

Voltage: ----- variable: 100 V~/ 115 V~/ 230 V~, 50–60 Hz
 Fuse, power supply: ----- 2 fuses, T 3.15 AL 250 V AC
 Power consumption: ----- 120 VA
 Applied part: ----- Type BF*
 Protection class: ----- Class I
 Dimensions (W x D x H): ----- 260 x 250 x 110 mm
 Net weight control unit: ----- 3.7 kg

Motor:

Motor coupling:----- Intra coupling ISO3964
 Motor speed: ----- 300 - 50,000 rpm**
 Max. Motor torque: ----- 6 Ncm
 Motor weight: ----- 0.280 kg
 Motor cable length----- 3 m

Pedal:

IP code (pedal) ----- IPX8

**Applied part of Type BF is the instrument used with the HighSurg 30.*

*** Optional available, Motor with a max. speed of 80,000 rpm and Motor 16 with integrated OTO-Drill Handpiece with 16'000 rpm.*

1.4 Ambient conditions

	Transport and storage:	Operation:
Relative humidity:	Max. 90%	Max. 80%
Temperature:	0°C – 60°C	10°C – 30°C
Atmospheric pressure:	700 hPa to 1060 hPa	800 hPa to 1060 hPa

1.5 Warranty coverage

Purchasing the HighSurg 30 entitles you to a 1-year warranty. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further **6 months**. Consumable parts are not covered by the warranty. Improper use or repair, or failure to observe these instructions, relieve us from any obligation arising from warranty provisions or other claims.

2 Explanation of symbols

	Important information		Autoclavable at 135°C
	Do not use if the packaging is damaged		Suitable for thermal disinfection
	Warning		Sterilized using ethylene oxide
	Manufacturer		Observe the instructions for use
	Type BF applied part Applied parts is the instrument		Symbol indicating the serial number with the date of manufacture (year/month).
	Do not reuse		Symbol indicating the order number.
	Biohazard		Symbol indicating the lot number.
	Socket for Motor 1		Socket for Motor 2
	Pedal		Containing phthalates
IPX8	Protection against continual submerging.		Certified by the Canadian Standards Association (CSA)
	Date of manufacture		Date of expiry
	CE symbol with notified body		Warning, heated up surfaces
	Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Valid local disposal regulations apply.		

EN

3 Safety information

Your safety, the safety of your team, and of course that of your patients is very important to us. It is therefore essential to bear the following information in mind:

Every use of the HighSurg 30 different to the product description defined in chapter “Intended use and operation”, causes risks for patients and trained personnel. If physical examinations and therapies are carried out without use of the devices then the devices must be removed from the place of treatment. Avoid any connection or close adjacency to other devices.

3.1 EMC Manufacturer's Declaration of Conformity

The use of (RF) Radio Frequency emitting devices and equipment as well as the occurrence of negative environmental factors in the close area of the HighSurg 30 may cause unexpected or adverse operation. Use only accessories and cables as specified in the product description. Further observe the EMC manufacturers declaration of conformity.

3.2 Integrated peristaltic tube pump

The integrated peristaltic pump is used to cool tissue in order to prevent damage to the tissue. It may only be operated with watery solutions, such as 0.9 % Sodium Chloride irrigation solution (Ref. 1707) or “Ringer” solution. Supplying medication using the integrated pump is expressly prohibited.

3.3 Modification and misuse

If modifications are made by the user/operator or a third party to the HighSurg 30 and the accessory equipment provided, or if these items are used by these individuals for a purpose other than that for which the items were designed (according to chapter 1.1 Intended use and operation), the manufacturer assumes no responsibility and the guarantee is void.

3.4 Essential requirements

The HighSurg 30 may only be operated by qualified and trained personnel.

Repairs may only be performed by authorized NOUVAG service technicians.

The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.

Prior to using the device, before startup, and before operation, the user must always ensure that the device and accessories are in good working order and are clean, sterile and operational.

Repairs may only be performed by authorized NOUVAG service technicians. Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims.

Use Nou-Clean spray for maintenance and care of the motor, handpieces and contra angles. Using other care products can result in malfunctions and/or cause the warranty to be revoked.

3.5 During use

The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Chapter 8 Cleaning, disinfection and sterilization).

At the choosing of the instrument the operator has to make sure it's biocompatible, according to EN ISO 10993.

Handpieces and contra angles may only be attached when the electronic motor is not running.

Never operate the clamping mechanism of the handpieces or contra angles while the system is running. This could result in instrument damage.

To prevent injury, never touch drill bits or burrs while they are still running.

Do not use device in the vicinity of flammable mixtures!

The user must ensure that as little friction heat as possible is generated during a procedure on a patient. High speeds and high application pressure may cause thermal necrosis of tissue.

4 Scope of delivery

Ref.	Description	Quantity
3360	HighSurg 30 control unit-----	1 unit
1510	Vario pedal; IPX8; electronic -----	1 unit
2099	Electronic motor 21, 3 m motor cable, max. speed 50,000 rpm-----	1 unit
 6024	Tube set, sterile, 3 m, single-use-----	1 unit
 6025	Tube set, sterile, 3 m, single-use (simultaneous coolant supply of two handpieces)-----	1 unit
1873	Clip set (10 pieces) for tube set attachment at motor cable-----	1 unit
1881	Clip set (3 pieces) for tube set attachment at handpieces and contra angles-----	1 unit
1707	Irrigation fluid; 0.9 % sodium chloride/water solution, 1 liter-----	1 unit
1770	Stand for irrigation fluid-----	1 unit
1170	Handpiece cradle-----	1 unit
1974	Spray nozzle for Nou-Clean spray; for care and maintenance of electronic motor -----	1 unit
1958	Spray nozzle for Nou-Clean spray; for care and maintenance of instruments -----	1 unit
31666	User Manual HighSurg 30 on CD-ROM -----	1 unit

In line with regulations pertaining to hazardous materials, the following items are not delivered with the control unit and must be ordered separately:

1984 Nou-Clean maintenance and care spray ----- 1 unit


Optional:

2098 Electronic motor 21 High Speed, 3 m motor cable, max. speed 80,000 rpm ----- 1 unit
 1909 OTO-Drill handpiece with integrated Motor 16, length of cable 3 m, max. 16'000 rpm ----- 1 unit

EN

5 Device overview


Front view

1. Pedal socket, device rear
2. Vario pedal
3. Handpiece (not included in delivery)
4. Clip for tubing set attachment at handpieces and contra angles
5. Handpiece cradle
6. Electronic motor (delivery includes 1 motor)
7. Clip for tubing set attachment at motor cable
8. 1 indicator light for each motor
9. 2 Motor sockets
10. Operating panel
11. Display
12. Release key for tubing set bracket
13. Peristaltic pump

14. Tubing set
15. Stand for irrigation fluid bottle
16. Drip chamber
17. Bleed valve
18. Bottle with irrigation fluid
19. Type plate with type designation, reference number, serial number, information on power supply and device fuse.
20. Potential equalization
21. Power plug socket
22. Main switch
23. Fuse compartment
24. Spray nozzle for maintenance of electronic motor
25. Spray nozzle for handpiece maintenance

Rear view

6 Startup

6.1 Device setup

EN

- Installation layout


- Place the HighSurg 30 and all required accessories and instruments on an even, non-slip surface and make sure you have good access to all controls.
- Do not allow the operating range of the device (including the cable and contra angle) to be compromised by limiting factors.
- The system display must be fully visible at all times.
- The pedal must be placed within stepping distance between the patient and the surgeon.
- It must be explicitly ensured that no objects can fall on the pedal.
- The power plug at the rear of the device must be accessible at all times.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

6.2 Connection to the power supply

Before plugging the power cable into the power socket for the first time, you must check the supply voltage setting next to the power switch.

If the voltage shown does not correspond to the local mains voltage, the grey fuse holder must be set to the correct voltage:

- Unplug the power cable.
- Use a screwdriver to open the fuse slot.
- Remove the fuse holder.
- Remove the grey fuse holder and reinsert it so that the local mains voltage setting is shown in the small window.
- Slide fuse holder back in and close the fuse slot.
- Check the mains voltage shown on the fuse slot.
- Plug the power cable back into the device.

In order to prevent the risk of electric shock, the device may only be connected to a power network with a PE protective ground conductor.

6.3 Device preparation

1. Sterilize the motor (the motor is not sterile on delivery). If the motor has already been sterilized: when removing the motor from the sterile packaging, ensure that the sterile packaging is not damaged and that the sterility indicator confirms sterility (if no sterility indicator is provided, the sterile packaging must at least show the date on which the shelf life of the sterile item is due to expire).
2. Insert the stand for the irrigation fluid into the stand holder.
3. Plug the motor plug of the electronic motor into one of the motor sockets.
4. Where appropriate, plug the motor plug of a second electronic motor into a motor socket.

- If you want to use a handpiece with the 80,000 rpm high speed motor and change its parameter, then the motor has to be plugged in one of the motor sockets. Otherwise there will be only displayed the values of the standard 50,000 rpm motor.
- The same is with the OTO-Drill handpiece with its integrated 16,000 rpm motor.

Note: Always plug in the motor first, then change the parameters of the selected handpiece.

5. Plug the pedal plug into the pedal socket at the rear of the control unit.
6. Attach the sterilized handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor until it clicks into position and make sure it is secure by moving it slightly in the opposite direction.
7. Assemble the tubing set: Decide whether you will use tubing set Ref. 6024, for the cooling of a single contra angle or if you need to use tubing set Ref. 6025 with an integrated 3 way cock for cooling two contra angles at the same time when two motors are used.

Use only Nouvag tube sets Ref. 6024 and 6025, otherwise the correct function cannot be guaranteed.


Check the expiry date of the tubing set and ensure that the packaging is not damaged. Using non-sterile tubing sets can result in serious infection and, in extreme cases, can be fatal.

When inserting the tubing set, watch the arrow on the cover of the pump compartment. It indicates the flow direction of the cooling liquid.

Do not regulate the amount of irrigation fluid using the roller clamp on the tubing set; with the HighSurg 30, this is regulated instead using the integrated pump. For this reason, make sure to open the roller clamp as far as it will go (please refer to 7.4.5 Setting the pump supply quantity).

- A) Press release key for tubing set bracket (on top of the control unit) to open the pump.
- B) The compartment with the integrated tubing bracket opens.
- C) Place the tubing set into the tubing bracket provided in such a way that the part of the tubing set with the spike exits the pump towards the rear of the device. Make sure the tubing is secure.
- D) With the tubing set inserted, press the compartment downwards until it clicks into place.


8. Insert the spike at the end of the tubing set into the irrigation fluid bottle and hang the bottle onto the stand.
9. Open the roller clamp on the tubing set as far as it will go.
10. Open the bleed valve beneath the drip chamber.
11. Connect the control unit to the power socket.

Ensure that the operating voltage setting corresponds to the local mains voltage.

6.4 Assembly of external irrigation system

- A) Secure motor cable clips and handpiece clips to the irrigation tube.
 - B) Connect the tube set with the external cooling tube of the Instrument (Example: Shaver blade).
 - C) Secure clips to the handpiece and to the motor cable.
- If required, secure additional clips to the motor cable.


7 Operation

7.1 Switching the device on and off

The power switch “I/O” (at the rear) is used to switch the control unit on and off.
The device can be switched off at any time irrespective of any procedure for switching off the device.


EN

7.2 Overview: Control elements on the operation panel

- A) **Display:** Shows the operating values (See chapter “7.3 Overview: Standard display”).
- B) **“Program” key:** Selection of program 1 to 10. There are multiple programs for each motor selectable. (To facilitate the selection of programs, the selectable programs can be limited to the number you need. Refer to “7.7 Configuration Menu, Level 1, Programs/Number of programs”).
- C) **“Motor 1/2” key:** Switching between the connected motors (the green indicator lights beside the motor sockets indicate the active motor). Pressing the Motor key longer changes direction of rotation.
- D) **“Handpiece” key:** Choice of handpieces or contra angle (can be individually deactivated, refer to “7.7 Configuration Menu”).
- E) **“Enter” key:** Used in the configuration menu (refer to “7.7 Configuration Menu”).
- F) **“Select” keys:**
 - By pressing the left “Select” key “**Select ▲**” the software version is shown.
 - By pressing the right key “**Select ▼**”, the type of connected motor is displayed.
 - By pressing both “Select” keys at the same time “**Select ▲ + ▼**” the programs will be reset to factory default settings.
 - The Select keys serve in the configuration menu to change values and parameters
“▲”: value adjustment (upwards) “▼”: value adjustment (downwards)
- G) **“Speed” keys:**
 - Restrict the maximum speed that can be selected using the pedal.
“▲”: increases the maximum speed “▼”: reduces the maximum speed
 - By pressing both “Speed” keys at the same time “**Speed ▲ + ▼**” the handpiece calibration will be started (refer to “7.4.2 Calibrating handpieces”).
- H) **“Torque/Frequency” keys:**
 - Restricts the maximum torque.
“▲”: increases the maximum torque “▼”: reduces the maximum torque
 - With the shaver handpiece attached the key pair is used to set the oscillation period, if oscillation mode is selected.
- I) **“Pump” keys:**
 - Changing the pump flow rate that can be supplied using the pedal.
“▲”: increases the maximum supply quantity “▼”: reduces the maximum supply quantity
 - By pressing both “Pump” keys at the same time “**Pump ▲ + ▼**” the pump will be put on call, pressing again will switch it off.

7.3 Overview: standard display

- A) Information Line**
Information and error messages are displayed here. Display is illuminated red when error messages are displayed.
- B) Program**
Shows selected number of program for the active motor.
- C) Pump**
The numerical value shows the pump flow rate in percent and the drop symbol together with the on/off indication shows if the pump is in stand-by mode or switched off.
- D) Rotational direction of the motor**
The arrow indicates the rotational direction set for the motor. The rotation direction can be changed by pressing the button “” on the pedal or by long pressing of the “Motor 1/2” key at the operating panel. In oscillating Shaver mode both arrows are displayed at the same time, which symbolizes the oscillation mode of the shaver.
- E) Current speed**
At stand still of the motor the set maximum speed is displayed. As the motor starts to run when the pedal is pressed, the evolving speed is shown in real-time.
- F) Date**
- G) Clock**
- H) Motor**
Shows selected motor → which is confirmed by the green LEDs next to each of the motor sockets.
- I) Name of the handpiece or corresponding transmission ratio**
Shows name of the handpiece used or the selected transmission ratio. (See also 7.4.1, “Selecting handpiece or transmission ratio”)
- J) Speed range**
Shows speed range of the handpiece used.
- K) Maximum torque**
Shows maximum torque setting.
- L) Current torque**
Bar graph providing a graphical representation of the current torque. All bars active means max. torque reached.

The pump does not begin to operate until the motor has been activated by pressing the pedal.


7.4 Adjusting the programs

Values for operation settings depend on the connected handpiece or contra angle as well as the task to be performed.

7.4.1 Step 1: Selecting handpiece or transmission ratio

Depending on the handpiece or contra angle attached to the motor, the corresponding transmission ratio has to be adapted accordingly by using the handpiece key.

Press the “Handpiece” key repeatedly until the name of the required handpiece with the corresponding transmission ratio is shown on the display. When the key is pressed constantly the handpieces or contra angles will be shown in fast forward mode.

Table of possible Handpieces, available for the 50,000 rpm motor

Name of Handpiece with transmission ratio	Handpiece name on Display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm
Handpiece, 1:5	1 : 5	1500	240,000	fix 1	
Handpiece, 1:3	1 : 3	900	150,000	1	2
Handpiece, 1:2	1 : 2	600	100,000	1	2
Handpiece, 1:1	1 : 1	300	50,000	1	6
Micro saws (Com., Osc., Sag.)	Micro Saw	fix 15,000		fix 6	
Dermatome	Dermatome	fix 14,000		fix 6	
Kirschner Handpiece	Kirschner	500	2800	fix 48	
Tattooing Handpiece	Tattoo Hp	fix (9000/10,000/11,000/12,000)		fix 6	
Jacobs Chuck	Jacobs Chuck	200	2600	fix 60	
Perforator	Perforator	80	900	20	80
Craniotome	Craniotome	1000	50,000	fix 6	
Handpiece 4:1	4:1	200	12,000	1	18

Please note: The Shaver parameters are listed in a separate table on the next page

Table of possible handpieces, available for the 80,000 rpm motor

Name of Handpiece with transmission ratio	Handpiece name on Display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm
Handpiece, 1:5	1 : 5	1500	240,000	fix 1	
Handpiece, 1:3	1 : 3	900	150,000	fix 1	
Handpiece, 1:2	1 : 2	600	100,000	fix 2	
Handpiece, 1:1	1 : 1	300	80,000	fix 3	
Jacobs Chuck	Jacobs Chuck	80	2100	fix 35	
Perforator	Perforator	80	1200	20	60
Craniotome	Craniotome	1000	60,000	fix 3	

Table of Handpiece OTO-Drill with integrated 16,000 rpm -motor

Motor with integrated Handpiece OTO-Drill	Handpiece name on Display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm
OTO-Drill	OTO Drill	1000	16,000	fix 1	

Shaver functions table, available for the 50,000 and 80,000 rpm motor

Name of handpiece with transmission ratio	Display	Speed min. rpm	Speed max. rpm	Oscillation period from ... seconds	Oscillation period to ... seconds
Continuous Shaver	Cont. Shaver	300	6000		
Oscillating Shaver	Osc. Shaver	300	5000	0.20	3.00

The oscillation period for oscillating shavers is determined by the speed. If the oscillation period is set too low for the set speed, it switches automatically to the next possible level.

This ensures that the shaver knife makes a full rotation for the speed setting.

Handpieces that don't belong to one's own assortment can be deactivated in the Configuration Menu. (refer to "7.7 Configuration Menu") That makes finding the right handpiece in the handpiece list easier, because the list of available handpieces becomes shorter.

7.4.2 Step 2: Calibrating handpieces

To make sure the displayed parameters correspond with the actual, measurable parameters of the handpiece or contra angle, it is recommended to calibrate each handpiece or contra angle on a regular basis.


It's a procedure as simple as it is important to guarantee safety and precision of each handpiece or contra angle being used.

After you take care of all prior preparations such as sterilization, maintenance and care of handpieces, device preparation and the selection of the handpiece of use, the calibration procedure can be performed.

The calibration of handpieces or contra angles guarantees accurate torque. Due to wear as well as varying lubrication of the handpieces or contra angles, and lack of maintenance and care, the distribution of torque can vary widely.

1. Select the correct handpiece, corresponding to the one mounted on the motor and double check on the display for correct selection.
2. Hold motor with mounted handpiece or contra angle in your hand at safe distance to your body.
3. Press both "Speed" keys at the same time (Speed ▲ + ▼). "Testing the handpiece XX" is displayed.

4. Motor and handpiece start running and pass several speed cycles up to maximum speed.
5. After a tone is emitted the calibration is finished. The display shows the message "Handpiece is OK".

If a handpiece does not work within the parameters established by calibration even after cleaning and lubrication, the device displays a red-illuminated error message, "Handpiece XX is faulty".

This indicates soiling, wear or a technical defect. Such handpieces must be repaired or replaced.

The torque response of all handpieces is checked during calibration.

With the 1:1 handpiece, control of the device is also adjusted to altered handpiece conditions to remain within the torque response tolerance.

7.4.3 Step 3: Setting the speed

The possible speed range depends on the attached handpiece. The maximum speed within this speed range can be restricted to the required value.

Using the pedal, the speed can be varied from the minimum speed up to the maximum speed as set.

Setting the speed:

Press the «Speed» keys “▲” to increase or “▼” to decrease the maximum speed. When key is pressed continuously the speeds will be shown in fast forward mode.

The following handpieces run only with one specific speed, it cannot be changed. The specific rpm values of all handpieces are shown in the table of chapter 7.4.1.

- Micro Saws (Oscillating Saw, Saggital Saw, Compass Saw)
- Dermatome
- Tattooing handpiece (4 steps)

7.4.4 Step 4: Setting the torque

Once the speed has been selected, the torque can be determined from the corresponding torque range.

Press the “Torque” keys “▲” to increase or “▼” to decrease the maximum torque. By pressing the key continuously the torques are shown in fast forward mode.

When the shaver handpiece is used in oscillation mode, the number of oscillations per second of the shaver knife (number of changes of direction of the rotating knife per second) is set with these two keys. In this function, the arrow displays both directions.

The following handpieces run on a fixed torque.

- Handpiece, 1:5 (1 Ncm)
- Micro Saws (6 Ncm)
- Dermatome (6 Ncm)
- Kirschner Handpiece (48 Ncm)
- Continuous Shaver (12/12 Ncm)*
- Oscillating Shaver (12/12 Ncm)*
- Tattooing Handpiece (6 Ncm)
- Jacobs Chuck (60/60 Ncm)*
- Craniotome (6/6 Ncm)*
- OTO-Drill with integrated motor (1 Ncm)

* The value before the slash refers to the use with the 50,000 rpm motor, the value after the slash to the 80,000 rpm motor.

7.4.5 Step 5: Setting the pump supply quantity

Press the “Pump” keys “▲” to increase or “▼” to decrease the pump supply quantity. By pressing the key continuously the values are fast forwarded.

Minimum and maximum of pump supply quantity as well as operating steps can be adjusted in the Configuration Menu (refer to «7.7 Configuration Menu, Parameter level 2, Pump»).

To activate or deactivate the pump, press both “Pump” keys at the same time, “Pump ▲ + ▼”, or use the foot switch .

7.5 Torque limiter

The device has an automatic torque limiter, which functions like a torque wrench.

If the attached instrument meets resistance, the torque increases as far as the maximum level, and the speed then reduces, to standstill if necessary. The torque on the instrument is maintained. If the load on the instrument reduces, the speed increases again up to the set maximum level.

This process can easily be seen on the display by means of the bar chart. The segments of the bar chart fill up as the resistance on the instrument increases. When the torque has reached the set maximum, that is, when all segments are visible, the speed reduces. As soon as the pressure on the instrument is reduced, the torque falls again. The bar chart on the display regresses and the instrument speed increases again.

7.6 Programme memory

At switch off of the HighSurg 30 all the settings made by the user are automatically stored in each program place. The display shows every setting of the selected program place. The following parameters are automatically stored:

- Handpiece/Transmission ratio
- Maximal speed
- Maximal torque
- Pump On/Off
- Pump performance
- Shaver oscillation period

To change a program go to the specific parameter and change it. When the device is switched off all parameters are automatically stored under the same program number.

In the configuration menu the number of storable programs can be limited for each motor separately.

In the configuration menu it can be chosen what the display shows at start-up of the device. Choose between Program 1 with Motor 1 or the last used setting before the device was switched off.

Number of available storing places for each motor and both motor connections

Motor-Choice	Motor connection 		Motor connection 	
	Number of available storing places	Factory default setting	Number of available storing places	Factory default setting
Motor 50,000 rpm	3 – 10	10	3 – 10	10
Motor 80,000 rpm	1 – 7	6	1 – 7	6
Motor 16,000 rpm, OTO-Drill	1 – 3	1	1 – 3	1

7.7 Configuration Menu

In the Configuration Menu the user can customize the device after his favor. The parameters are organized in two levels. "Parameter Level 1" comprises all basic information and possible settings that affect the device itself. Parameter Level 2 includes all information and settings that involve the surgeon's work. In the configuration menu the following information can be obtained or parameters can be set according to one's own needs:

- Software version
- Mainboard serial number
- Date format US ?
- Date
- Time
- Battery voltage
- Language, DE/EN
- Display brightness
- Number of programs for 50,000 rpm motor
- Number of programs for 80,000 rpm motor
- Number of programs for Motor 16 OTO-Drill
- Behavior after switching on
- Operating hours HighSurg 30
- Operating hours Motor 1
- Operating hours Motor 2
- Operating hours of irrigation pump
- Error-messages (the last 8)
- Activation of available handpieces (14 pos.)
- Enter Password
- Speed limitation for each handpiece
- Pump behavior (10 pos.)
- Motor behavior for 50,000 rpm motor
- Motor behavior for 80,000 rpm motor
- Motor behavior for Motor 16 OTO-Drill
- Reset to factory default

Be cautious when changing parameters. Unusual behavior of instruments while operating may provoke false reactions and jeopardize the patient. Every setting and the new behavior of the instrument has to be verified and tested.

Access to Configuration Menu:

- Press **"Enter"** key for 3 seconds until a long tone is emitted. The right arrow **"➤"** in the information line on the display indicates that you're in the Configuration Menu.

Enter

```
>Software
>Version          U1.13
PROGRAM 03 MOTOR 2
```

Changing certain parameters:

- Choose the desired parameter by pressing **"Select ▲"** or **"Select ▼"**.
- Press **"Enter"** to activate setting mode. The value of the parameter is now displayed in parentheses.
- To change the value press **"Select ▲"** or **"Select ▼"**.
- To confirm the changes made, press **"Enter"** for one sec., until a short tone is emitted.
- To abort the settings made, press **"Enter"** briefly, the setting returns to its previous value, without emitting a tone.

```
>Backlight
>brightness      (0..10) [02]
```

```
>Backlight
>brightness      (0..10) 10
```

Exit from Configuration Menu:

- To leave the Configuration Menu press **"Enter"** for 3 seconds, until a long tone is emitted.

Parameter Level 1

Group/Parameter	Rights	Factory	Definition
Software/Version	read	VX.XX	Shows current software version
Hardware/Serial number MB	read	92XXXXXX	Shows serial number of the main board
Date-Time/Date format US	read/write	no	Set to US Date format
Date-Time/Date	read/write	–	Change current date
Date-Time/Time	read/write	–	Change current time
Battery/Voltage	read	e.g. 3120 mV	State of charge of battery (Capacity about 3100 mV, Low 1800 mV)
Language/o = DE, 1 = EN	read/write	1	Change display language between German and English
Backlight/brightness (0..10)	read/write	9	Display brightness, variable: 0, ..., 10
Program/Number of programs for 50,000 rpm motor	read/write	10	Number of activated Programs: 3, ..., 10
Program/Number of programs for 80,000 rpm motor	read/write	6	Number of activated Programs: 1, ..., 7
Program/Number of programs for 16,000 rpm motor	read/write	1	Number of activated Programs: 1, ..., 3
Program/Power on at last program	Read/write	Yes	No: Display of Motor 1 with Program 1 Yes: Display of the last used program/motor
Operating hours/Control Unit	read	XX:XX:XX	Shows operating hours of HighSurg 30
Operating hours/Motor 1	read	XX:XX:XX	Shows operating hours of motor 1
Operating hours/Motor 2	read	XX:XX:XX	Shows operating hours of motor 2
Operating hours/Pump	read	XX:XX:XX	Shows operating hours of pump
Error memory/ 1 – 8	read	X	8 Error messages in chronological order. <i>For the list of error codes please refer to chapter 10.</i>

Handpiece activation	Name of handpiece on display	Choice	Factory default	Definition
Handpiece existing/HP 01	1 : 5	yes/no	yes	Deselect handpieces that do not belong to your product range by switching to "no". This will later shorten the scrolling list for "Handpieces". Otherwise you will later have to scroll down the complete list of available handpieces, each time you select a handpiece. <i>Because HP 04, 1:1 is permanently active, it's not listed in this table.</i>
Handpiece existing/HP 02	1 : 3	yes/no	yes	
Handpiece existing/HP 03	1 : 2	yes/no	yes	
Handpiece existing/HP 05	Com.Saw	yes/no	yes	
Handpiece existing/HP 06	Osc.Saw	yes/no	yes	
Handpiece existing/HP 07	Sag.Saw	yes/no	yes	
Handpiece existing/HP 08	Dermatome	yes/no	yes	
Handpiece existing/HP 09	Kirschner Hp	yes/no	yes	
Handpiece existing/HP 10	Tattoo Hp	yes/no	yes	
Handpiece existing/HP 11	Cont.Shaver	yes/no	yes	
Handpiece existing/HP 12	Osc.Shaver	yes/no	yes	
Handpiece existing/HP 13	Jacobs Chuck	yes/no	no	
Handpiece existing/HP 14	Perforator	yes/no	no	
Handpiece existing/HP 15	Craniotome	Yes/no	no	

Parameter Level 2

Values in Parameter Level 2 can only be adjusted after the password “9403” was entered. The password cannot be changed.
 Entering password: Press “Enter”, press “Select ▲ or ▼”, (for fast forward or backward keep key pressed).
 To confirm press “Enter” for 1 second, until a tone is emitted.

Handpiece Max. Speed	Name of handpiece on the display	Speed range rpm	Factory default	Definition
Handpiece max speed/HP 01	1 : 5	1500 – 240,000	240,000	Limit the maximum speed of your handpieces according to your own experience.
Handpiece max speed/HP 02	1 : 3	900 – 150,000	150,000	
Handpiece max speed/HP 03	1 : 2	600 – 100,000	100,000	
Handpiece max speed/HP 04	1 : 1	300 – 80,000	80,000	
Handpiece max speed/HP 07	Kirschner Hp	500 – 2800	2800	
Handpiece max speed/HP 08	Tattoo Hp	9000 – 12,000*	12,000	
Handpiece max speed/HP 09	Cont. Shaver	300 – 6000	6000	
Handpiece max speed/HP 10	Osc.Shaver	300 – 5000	5000	
Handpiece max speed/HP 11	Jacobs Chuck	200 – 2600	2600	
Handpiece max speed/HP 12	Perforator	80 – 1200	1200	
Handpiece max speed/HP 13	Craniotome	1000 – 60,000	60,000	
Handpiece max speed/HP 14	4 : 1	75 – 20,000	20,000	
Handpiece max speed/HP 15	OTO-Drill Hp	300 – 16,000	16,000	

* Speed selection in 4 steps: 9000, 10'000, 11'000, 12'000 rpm

Pump parameters	Range	Factory default	Definition
Pump/Backwards turn mode variable	No/Yes	Yes	The pressures in the tube set vary according to the pump speed. In “variable mode” this is considered, to prevent of spilling when pump is switched off in backwards mode.
Pump/Way backwards	1 – 100%	25%	Specify how far the pump turns backwards
Pump/Speed backwards	10 – 50%	33%	Specify how fast the pump has to turn backwards to prevent of spilling after the switch off of the motor.
Pump/Range 1 increment	1 – 10%	1%	Adjustment steps in section 1.
Pump/Range 1 end	5 – 50%	20%	Set the range where section 1 is active
Pump/Range 2 increment	1 – 10%	5%	Adjustment steps in section 2
Pump/Range 2 end	10 – 90%	50%	Set the range where section 2 is active.
Pump/Range 3 increment	1 – 10%	10%	Adjustment steps in section 3
Pump/Range 3 end	20 – 100%	100%	Set the range where section 3 is active.
Pump switch mode on/off	No/Yes	No	No: Pump starts simultaneously with the handpiece drive. Yes: Pump starts separately by pressing the foot switch. 

The HighSurg 30 can recognize the type of a plugged-in motor. This enables the adaption of the motors and their safe operation.

Motor 21, 50,000 rpm	Range	Factory default	Definition
Motor 21, 50,000 rpm/Min. speed	300 – 5000 rpm	300 rpm	Set the min. speed
Motor 21, 50,000 rpm /Max. speed	5000 – 50,000 rpm	50,000 rpm	Set max. speed of
Motor 21, 50,000 rpm/Ramp start	1 – 1000 ms/10,000 rpm	250 ms	Set acceleration time to 10,000 rpm
Motor 21, 50,000 rpm/Ramp stop	1 – 1000 ms/10,000 rpm	50 ms	Set breaking time from 10,000 – 0 rpm

Motor 21, 80,000 rpm	Range	Factory default	Definition
Motor 21, 80,000 rpm/Min. speed	300 – 5000 rpm	300 rpm	Set the min. speed of
Motor 21, 80,000 rpm /Max. speed	5000 – 80,000 rpm	80,000 rpm	Set max. speed of
Motor 21, 80,000 rpm/Ramp start	1 – 1000 ms/10,000 rpm	250 ms	Set acceleration time to 10,000 rpm
Motor 21, 80,000 rpm/Ramp stop	1 – 1000 ms/10,000 rpm	50 ms	Set breaking time from 10,000 – 0 rpm
Max. speed warning delay	0 – 1000 ms	1000 ms	0: Warning is switched off 1 – 1000: Delay time of warning when exceeding 50,000 rpm.


Motor 16, OTO-Drill	Range	Factory default	Definition
Motor 16, OTO-Drill/Min. speed	1000 – 5000 rpm	1000 rpm	Setting of the min. speed
Motor 16, OTO-Drill/Max. speed	5000 – 16,000 rpm	16,000 rpm	Setting of max. speed
Motor 16, OTO-Drill/Ramp start	1 – 1000 ms/10,000 rpm	250 ms	Set acceleration time to 10,000 rpm
Motor 16, OTO-Drill/Ramp stop	1 – 1000 ms/10,000 rpm	50 ms	Set breaking time from 10,000 – 0 rpm

Resetting to factory default	Choice	Factory default	Definition
Default value/Set default value	Yes/No	–	Resetting all parameters of the Configuration Menu to factory default.

- **Attention:** Upon resetting to factory default all parameters will appear with factory default values (except date and time and the operating hours counter).
- **Note:** For resetting the “programs” to factory default press both “Select” keys at the same time (“Select ▲ + ▼”). For this procedure you have to be out of the Configuration Menu.

7.8 Operation using the Vario pedal

- 1. **Carrying bracket**
The carrying bracket can be operated by foot (collapsible).
- 2. **Key **:
Pressing the key briefly: switches the pump on or off (check display).
Pressing the key longer: increases the pump speed (check display).

In the configuration menu, parameter Level 2 (Pump switch mode on/off) the pump can be set to work permanently without pressing the treadplate. Therefore only a press on the footswitch  is needed and the pump runs with the presetted value. Pressing the footswitch again stops the pump. In this case the pressing of the pedal has no effect.

- 3. **Key **:
Pressing the key briefly: switches the program (+1) (check display).
Pressing the key longer: switches the program (-1) (check display).
- 4. **Key **:
Pressing the key briefly: switches the rotational direction (check display).
Pressing the key longer: switches the motor (check display and indicator lights beside the motor sockets on the control unit).
- 5. **Treadplate**
With the treadplate pedal the motor speed is variably adjusted and the pump is activated.

Treadplate...	Motor:	Pump:	Pump permanently on: <i>If in the configuration menu, parameter Level 2 "Pump switch mode on/off" is activated and the pump is started by pressing the foot switch. </i>
... not pressed	Motor off	Pump off	Pump on, if pump "On" displayed (speed as set on the control unit)
... pressed gently	Motor runs slowly	Pump on, if pump "On" displayed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed as set on the control unit)
... pressed all the way down	Motor runs at maximum speed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed as set on the control unit)

For safety reasons, the unit can only be operated by pedal.

The speed of the following handpieces is fixed and cannot be changed.

- Micro saws
- Dermatome
- Tattooing Handpiece (4 steps)

7.9 Functional check

Prior to HighSurg 30 startup or use of accessory equipment, the user must always ensure that each individual component is in good working order, free from defects, and is clean, sterile and operational. All inscriptions on the device and its accessories must be readable and there must be no loose parts in the device. Once the device is switched on, depending on configuration the most recent settings entered appear on the display and the green LED for motor 1 lights up.

7.9.1 Electronic motor

The performance of the electronic motor is checked without a handpiece or contra angle. Nevertheless, the **1:1 handpiece must be activated via the «Handpiece» key** to check the maximum speed. Use the “Speed” selection keys to set the speed of the electronic motor to 50,000 rpm. Press the pedal treadplate; the electronic motor starts and accelerates to up to 50,000 rpm. When the treadplate is released, the electronic motor slows down again.

For motor types with a permitted speed of 80,000 rpm, the maximum speed must be set accordingly to 80,000 rpm and for the Motor 16 with integrated Handpiece to 16,000 rpm.

- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

7.9.2 Pump

Press the key  on the pedal briefly; the peristaltic pump is switched on, which is shown on the display by the symbol of a drop. Press the pedal treadplate; the peristaltic pump and the electronic motor start up. Water sprays from the irrigation needle on the contra angle.

In the configuration menu, parameter Level 2 (Pump switch mode on/off) the pump can be set to work permanently without pressing the treadplate. Therefore only a press on the footswitch  is needed and the pump runs with the presetted value. Pressing the footswitch again stops the pump. In this case the pressing of the pedal has no effect.

7.9.3 Rotational direction of the electronic motor

Press the key  on the pedal briefly; the rotational direction of the electronic motor changes. Press the pedal treadplate; the electronic motor rotates to the left and a continuous tone is emitted. Release the treadplate; the electronic motor ceases to operate and the tone is no longer heard. By pressing the motor key again, the rotational direction is switched back to right rotation, which is shown on the display by the symbol of a changing arrow.

7.9.4 Program

The required program can be set by repeatedly pressing the key  on the pedal.

8 Cleaning, disinfection and sterilization

The instructions described here are intended for the parts supplied in the set. The cleaning, disinfection and sterilization instructions for extensions and accessories are described in their respective operating instructions.

EN

The following points in particular are important with regard to caring for the material:

- Perform cleaning, disinfection and sterilization after every treatment.
- Always autoclave the material in sterilization packaging.
- Make sure that sterilization packaging is no more than 80 % full.
- Always autoclave the material at 135°C for at least 5 minutes.
- If sterilized material is not used immediately, the material packaging must be labeled with the sterilization date.
- Nouvag AG recommends including a sterility indicator.

8.1 Control unit and pedal



The control unit and pedal do not come into contact with the patient.
Wipe the outside using micro-biologically tested surface disinfectant or a 70 % isopropyl alcohol solution.
The front plate of the control unit is sealed for this purpose and can be wiped clean.

8.2 OTO-Drill handpiece with integrated Electronic Motor 16

The OTO-Drill handpiece is can be ordered optionally. The reprocessing instructions are contained in the operation manual, delivered with the product.

8.3 Electronic motor 21

Reprocessing restrictions	Frequent reprocessing has only a limited impact on the electronic motor. The end of the product service life is normally determined by wear and damage through use. The Electronic Motor 21 and the Electronic Motor 21 High Speed are designed for 250 sterilization cycles.
INSTRUCTIONS	
At the location of use	Remove surface soiling with a disposable cloth/paper towel.
Storage and transport	No special requirements. Due to the risk of drying and corrosion, reprocessing must be performed without undue delay.
Preparation for cleaning	Remove soiling from the electronic motor with a disposable cloth/paper towel. Unscrew the motor cap, remove the cable and unscrew the handpiece carrier. 
Automatic cleaning and disinfection	Equipment: Washer-disinfector with a special load carrier that ensures motors are connected to the washer-disinfector and channels are rinsed. Rinse the motor from the front. Only use a neutral cleaning agent for this purpose. Place the electronic motor in the load carrier (ensure that the channels can be rinsed). Place the motor cap and cable, and handpiece holder, in the basket. Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with fully deionized water. Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection. When removing the electronic motor, check the motor cap and cable, and handpiece holder, to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.
Manual cleaning	Equipment: Neutral cleaning agent, soft brush, running demineralized water (< 38°C) Procedure: Rinse off and brush away surface soiling on the electronic motor, motor cap and cable, and handpiece holder. Use a brush to apply cleaning agent to all surfaces and gaps. Rinse the electronic motor, motor cap and cable, and handpiece holder thoroughly under running water.
Manual disinfection	For manual disinfection, wipe the electronic motor, the motor cap and the plug and cable underneath it, and the handpiece holder, with a certified non-chlorine disinfectant.
Drying	If a drying program is not provided by the washer-disinfector, the electronic motor must be dried manually or in a drying cabinet. The handpiece holder must then be screwed back onto the motor.

Inspection and maintenance  Ref. 1974	Perform a visual inspection to check for damage, corrosion and wear. Connect Motor with the spray nozzle (Ref. 1974) by screwing it on and spray motor for about 3 seconds and briefly wipe with a moist cloth (see instructions on spray can). Once thorough spraying has been completed, screw the cable and motor cap back onto the electronic motor. 
Packaging	Individual: Pack the electronic motor in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. Nouvag AG recommends including a sterility indicator. Sets: Sort electronic motors on trays intended for this purpose or place them on all-purpose sterilization trays.
Sterilization	The steam sterilization is carried out according to sterilization cycle B (EN 13060 fractionated pre-vacuum)." Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the electronic motor to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards. * Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.
Storage	If the sterilized electronic motor is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Including a sterility indicator is recommended.

The effectiveness of the sterilization instructions provided above for electronic motor 21 and the handpiece has been validated by Nouvag AG. The user is responsible for ensuring that the sterilization procedure performed achieves the required result. This requires validation and routine monitoring of the procedure. The staff member who completes the procedure bears sole responsibility for any deviation on his part from the instructions provided. Deviations necessitate revalidation of the effectiveness of the procedure as well as of the technical resilience of the motor with regard to the modified sterilization process.

8.4 Tubing set, Ref. 6024 and Ref. 6025

- Single-use tubing sets 6024 and 6025 may not be reused.
- Tubing sets must be disposed of properly after use!
- Do not use tubing sets when packaging is already open or damaged!
- Do not use tubing sets when expiry date has passed.
- Use only Nouvag tubing sets with Ref. 6024 and Ref. 6025!

Sterility cannot be guaranteed by reusing and re-sterilization of tubing sets. The characteristics of the material change in a manner that can result in failure of the system. This may result in serious infections or even patient death worst-case.

8.5 Handpiece cradle


Soiled handpiece cradles are cleaned using a neutral cleaning agent and then sterilized in accordance with the same instructions as for electronic motor 21.

9 Maintenance

9.1 Replacing the control unit fuse

Users can replace faulty control unit fuses themselves. These are located at the rear of the device in the fuse slot beside the power switch:

- Unplug the power plug.
- Open the fuse slot using a screw driver.
- Replace the faulty fuse T 3.15 AL 250 V AC.
- Slide the fuse holder back in and close the fuse slot.
- Check the mains voltage shown on the fuse slot.
- Plug in the power plug again.

1. Fuse slot locking mechanism
2. Display window for voltage setting
3. Fuse slot
4. Fuse 1
5. Fuse 2

9.2 Safety inspections

The essential requirements have been defined and within the risk analysis assessed. The approved results have been filed in the risk management deposited with the manufacturer.

The performance of safety inspections on medical devices is required by law in several countries. The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

The STI (Safety technical inspection) for the HighSurg 30 shall be executed every 2 years by authorized experts. Results shall be documented.

NOUVAG AG offers a safety inspection service for its customers. Addresses can be found in the appendix of this operation manual under "Service centers". For further information please contact our technical service department.

Further international service centers are listed on the Nouvag website:

www.nouvag.com > Service > Service centers

10 Malfunctions and troubleshooting

EN

Malfunction	Cause	Solution	Refer to operating instructions
Device is not operational	Control unit is not switched on	Set the power switch "I/O" to "I"	7.1 Switching the device on and off
	Power connection not established	Connect the control unit to the mains power supply	6.2 Connection to the power supply
	Incorrect operating voltage	Check the mains voltage	6.2 Connection to the power supply
	Faulty fuse	Replace the fuse	9.1 Replacing the control unit fuse
Motor does not run	Motor not switched on	Switch on the motor using the treadplate	7.8 Operation using the Vario pedal
	Incorrect motor active	Switching to the other motor, using the Vario pedal	7.8 Operation using the Vario pedal
	Motor not connected	Connect the motor cable to the control unit	5.0 Device overview 6.2 Connection to the power supply
	Handpiece or contra angle not correctly assembled	Press the handpiece firmly onto the electronic motor until it clicks into position and check that it is secure by moving it slightly in the opposite direction.	6.3 Device preparation
No irrigation fluid for instrument	Peristaltic pump not switched on	Switch on the peristaltic pump	7.8 Operation using the Vario pedal
	Tubing set incorrectly inserted	Insert tubing set correctly (note the direction)	6.3 Device preparation
	Tubing set clogged/crusted matter visible	Replace the tubing set	6.3 Device preparation
	Bottle with sodium chloride solution not ventilated	Open the ventilation filter at the drip chamber	6.3 Device preparation
	Tubing set is dripping	Replace the tubing set	6.3 Device preparation
	Roller clamp of tubing set is closed	Open roller clamp all the way	6.3 Device preparation
Pedal is not operational	Pedal not connected	Connect the pedal to the control unit	6.3 Device preparation
	Incorrect operation	Check operating instructions	7.8 Operation using the Vario pedal

If a fault cannot be rectified, please contact your supplier or an authorized service centre. The addresses are provided on the last page of the operating instructions.

HighSurg 30, Error-messages on Display		
Error-messages/Error Code	Cause	Solution
Basic Initialization/ W00	First Initialization	
Set default value/ W01	Parameter reset to default value	
Memory error/ E02	System Error	Send Control Unit to Service Center.
Handling error/ E03	System Error	Send Control Unit to Service Center.
Program SW error/ E04	System Error	Send Control Unit to Service Center.
UserConfig SW error/ E05	System Error	Send Control Unit to Service Center.
Display error/ E06	System Error	Send Control Unit to Service Center.
Pump error/ E07	System Error	Send Control Unit to Service Center.
Storing factory settings/ User Config & Program	Message while default values of parameters and programs are saved to the NOU-dongle.	
Storing factory settings/ Program	Message while default values of programs are stored.	
Pedal not connected/ E10	a) Pedal is not plugged in. b) Plug or cable is defective.	a) Plug in Pedal b) Send Control Unit and pedal to Service Center.
Pedal test mode/ W11	Pedal test mode switched on	Switch off device for 5 seconds, than switch on again.
Keyboard test mode/ W12	Keyboard test mode switched on	Switch off device for 5 seconds, than switch on again.
No motor connected/ E13	a) No motor connected b) Motor, motor cable, motor plug or Control Unit are is defective	a) Plug in motor b) Send motor and Control Unit to Service Center.
Motor 2 not connected/ E14	a) Motor 2 is selected but no motor plugged in b) Motor connected to motor socket 2, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in motor b) Send motor and Control Unit to Service Center.
Motor 1 not connected/ E15	a) Motor 1 is selected but no motor plugged in. b) Motor connected to motor socket 1, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in motor b) Send motor and Control Unit to Service Center.
Unknown motor 2/ E16	a) Motor 2 is selected but wrong motor is plugged in. b) Right motor is connected to motor socket 2, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in correct motor b) Send motor and Control Unit to Service Center.
Unknown motor 1/ E17	a) Motor 2 is selected but wrong motor is plugged in.	a) Plug in correct motor. b) Send motor and Control Unit to Service

	b) Right motor is connected to motor socket 1, but motor, motor cable, motor plug or Control Unit is defective.	Center.
Hp not allowed for motor 2/ E18	The Handpiece may not be operated with the motor plugged in at the motor socket 2.	Choose appropriate Handpiece or change to correct motor type.
Hp not allowed for motor 1/ E19	The Handpiece may not be operated with the motor plugged in at the motor socket 1.	Choose appropriate Handpiece or change to correct motor type.
Pump is open/ E20	Motor is not working when pump compartment is open to prevent of injuries.	Close pump compartment.
Motor or pump test mode/ W21	Motor or pump test mode is switched on.	Switch off device for 5 seconds, than switch on again.
Pedal locked/ W26, pedal let go	If pedal is pressed at switch on procedure, pedal will not work.	Release pedal for one second.
Battery is almost empty/ W27, continue with „Enter“	Battery is almost empty.	- After pressing „Enter“, work can be continued. - Send Control Unit as soon as possible to Service Center.
Watch Error XX/ E28, continue with „Enter“	a) Clock at Control Unit is defective. b) Device was switched on after battery change but clock is not set yet.	a) After pressing „Enter“, work can be continued but Control Unit has to be sent to Service Center as soon as possible. b) Set date and time.
Handpiece XX is faulty/ E29	At calibrating or testing handpiece was overexposed to high torque.	- Clean handpiece and spray it thoroughly with Nou-Clean spray. - If message is still displayed after test procedure, handpiece has to be sent to Service Center.
Handpiece XX is Ok!	Tested handpiece is OK.	
Testing the handpiece XX	Handpiece is testing.	
Check max. speed handpiece/ W30, wait 1 second	Message is displayed when motor with 80,000 rpm is selected and max. speed is set to 50,000 rpm. Check Handpiece for max speed clearance before you proceed.	Check handpiece and wait for the message to disappear. Motor starts if the pedal still is pressed.
NOU-Dongle is plugged in	Message is displayed for 1 second after NOU-dongle was plugged in.	

The red highlighted error messages are also illuminated red on the control device display. The other messages are for information and do not require any action by the user.

11 Spare parts list with order numbers

Accessories	Ref.
Clip set large CL, for attachment to the handpiece, package with 3 units-----	1881
Clip set motor cable, for attachment to the motor cable, package with 10 units-----	1873
Single-use tubing set, 3 m, sterile, 10-unit pack-----	6024
Single-use tubing set with integrated 3 way cock, 3 m, sterile, 10-unit pack-----	6025
Irrigation fluid, 0.9 % sodium chloride, 1 liter-----	1707
Nou-Clean clean and care spray-----	1984
Spray nozzle attachment for surgical instrument e-coupling-----	1958
Spray nozzle attachment for electronic motor 21-----	1974

To order additional parts, please contact our customer service department.

User Manual for HighSurg 30----- 31666

No instructions for use in paper form are enclosed with this product. Instructions for use in PDF format are enclosed on CD-ROM and require a CD-ROM drive and Adobe Acrobat Reader software to enable the instructions for use to be displayed or printed out. If you prefer the instructions for use in paper form, it must be requested via the address or website stated on the type plate.

12 Information on disposal

When disposing of the device, device parts and accessories, the regulations prescribed by law must be observed.

Do not dispose of devices with household waste!

To ensure environmental protection, old devices can be returned to the dealer or manufacturer.

Motors that have reached the end of their service life may not be disposed of with household waste. Motors must be sterilized before disposal. Please observe currently valid national disposal regulations for infectious waste.

Contaminated single-use tubing sets are subject to specific disposal requirements. Please observe currently valid national disposal regulations for infectious waste.

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
Electromagnetic compatibility (EMC)

Remark:

The **Product** subsequently referred to herein always denotes the **HighSurg 30**.

Changes or modifications to this product not expressly approved by the manufacturer may result in increased emissions or decreased immunity performance of the product and could cause EMC issues with this or other equipment. This product is designed and tested to comply with applicable regulations regarding EMC and shall be installed and put into service according to the EMC information stated as follows.

WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the product may cause unexpected or adverse operation.

WARNING

The product shall not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the product shall be tested to verify normal operation in the configuration in which it is being used.

Compliant Cables and Accessories

WARNING

The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the product.

The table below lists cables, transducers, and other applicable accessories for which the manufacturer claims EMC compliance.

NOTE: Any supplied accessories that do not affect EMC compliance are not listed.

Description	Length max.
Electronicmotor 21	2.9m
Vario Footpedal IP68	2.9m

Guidance and manufacturer's declaration – electromagnetic emissions		
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Product requires continued operation during power mains interruptions, it is recommended that the Product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity for not life support equipment			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V rms 150 kHz to 80 MHz outside ISM bands	10 V rms 150 kHz to 80 MHz outside ISM bands	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
			where P is the maximum output power rating in the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a	Fixed strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Product is used exceeds the applicable RF compliance level above, the Product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Product.		
b	over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended separation distances between portable and mobile RF communications equipment and the not life support equipment			
The Product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Product as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the higher frequency range applies.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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For all Nuvag servicecenters worldwide please check:

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Nouvag AG • St.Gallerstr. 23–25 • CH-9403 Goldach
Nouvag GmbH • Schulthaißstrasse 15 • D-78462 Konstanz